

Dermatology

[Liquid Nitrogen Cryosurgery Instruments](#)

[Tip Attachments/Miscellaneous](#)

[Electrosurgical Unit](#)

[Smoke Evacuator](#)

[Bipolar Accessories](#)

[Disposable Biopsy Punch](#)

[Dermabraders](#)

[Instruments](#)

Cryosurgery

Nitrospray Plus

Trigger Operated, Closed System,
Liquid Nitrogen Cryosurgical Unit

- Improved Seal⁶
Eliminates o-ring positioning.
- Grooved Head
Easy to use, comfortable to hold.
- Curved Trigger
Unobstructed view of treatment site.
- Notched Head Grip
Easy to remove/install cap.
- Three (3) year warranty

30-Day
Unconditional
Warranty

Static holding time:
8 - 12 hours

Intermittent use:
3 - 6 hours of patient treatment time



1006060 Nitrospray Plus
(16oz fill capacity)
Includes 5 spray tips,
tip protector and
base holder



1006065 Nitrospray Plus Lite
(10oz fill capacity)
Includes 5 spray tips,
tip protector and
base holder



6. Pressure seal ring retrofit kit (1006066) to eliminate o-ring on units prior to RP-11061 and LP-10952 is available. Not to be used with RP-10454 to RP-10561. Contact Customer Service for assistance.

Cryosurgery - Tip Attachments

Spray Tips

(shown actual size)



Spray Tip, 16 gauge • 1006516
(purple)



Spray Tip, 17 gauge • 1006517
(white)



Spray Tip, 18 gauge • 1006518
(pink)



Spray Tip, 19 gauge • 1006519
(brown)



Spray Tip, 20 gauge • 1006520
(yellow)

Spray Tips, Set of Five • 1006525
one of each size (16, 17, 18, 19, 20 gauge)

Cutaneous Probes*

(not actual size)



Cutaneous Probe, 2mm • 1006602



Cutaneous Probe, 4mm • 1006604



Cutaneous Probe, 6mm • 1006606



Cutaneous Probe, 8mm • 1006608



Cutaneous Probe, 10mm • 1006610



Cutaneous Probe, Sharp • 1006620

Cutaneous Probes, Set of Five • 1006615
(2mm, 4mm, 6mm, 8mm, 10mm)

*Reusable vent tube provided with each probe.

Cryocones

(not actual size)



Cryocone, 2.5mm • 1006625



Cryocone, 5mm • 1006635



Cryocone, 10mm • 1006640



Cryocone, 20mm • 1006650




Cryocone, 25mm • 1006665

Cryocones, Set of Five • 1006670
(2.5mm, 5mm, 10mm, 20mm, 25mm)

Miscellaneous

(not actual size)



Bent Spray Extension • 1006530



Spray Tip Protector • 1006527



Needle Extender, Rigid • 1006700



Brymill Accessory Adapter • 1006521

1006730 Replacement Bulb/Hose for Original Nitrospray

1006735 2" Needle Extender



**Model 2001e
Electrosurgical Unit**

Full-Featured Electrosurgery
for the Office

Designed with the alternate
site in mind!

Increased capabilities beyond
coagulation-only units!

9006501 Model 2001e Electrosurgical Unit
(110/240VAC)
(Includes: hospital-grade
power cords (110V and 220V),
disposable hand-control
pencil, disposable patient pad,
clamp cord assembly and manual)

See page 2 for a full description, specifications
and warranty information.



**Evacuator/Air Filtration
System**

- 9006300 Premier Evacuator (110V)
Includes One (1) Filter

- 9006298 Premier Evacuator (220V)

- 9006290 Replacement Filter Packs (5 ea/case)

- 9006292 Reusable Smoke Tube
(for accessible sites)
Length: 8 7/8", Diameter: 3/4"

- 9006296 Evacuator Footswitch (not shown)

- 9006297 Hands-Free Tube Holder (not shown)
Eliminates the need of an assistant,
while providing a full range of motion.

See page 3 for a full description, specifications and
warranty information.

Electrosurgery

Bipolar Coagulating Forceps

Further expand the capability of the 2001e Electrosurgical Unit with Premier Bipolar Forceps. Fully autoclavable, with a standard two-prong connector, these forceps are compatible with other generators and bipolar cord sets.

Must be used with foot control (refer to 9006525 Dual Foot Controller below).



Because the output signal is contained between the forceps tips, no grounding pad is required to complete the electrical circuit.

This makes bipolar the preferred mode of coagulation when doing simple procedures around the head, face and upper chest.

As with any electrosurgical procedure, extreme caution should be used on patients with cardiac pacemakers, diagnostic monitoring systems or infusion pumps.

(Illustration above for descriptive purposes only. Not representative of actual procedure.)



9006308 5 1/2" Straight



9006310 7 1/2" Bayonet, Fine Tip



9006309 5" Curved



9006314 7" Straight, Fine Tip

Contact Premier Customer Service for the availability of additional forcep patterns

9006313 Standard Bipolar Adapter Cord, reusable (Not Shown)
Banana plug male connectors are compatible with most manufacturers' generators. Accepts forceps with male 2 pin connectors.

Dual Foot Controller


Fluid resistant, dual operating mode conveniently selects either cut or coagulation operating modes. Standard color-coded pedals included.



9006525 Dual Foot Controller


(For use with the Premier 2001e generator only.)

Disposable Biopsy Punch



Uni-Punch®

Sterile, Seamless, Razor Sharp Blades



Box of 25

- 9033509 1.5mm
- 9033502 2.0mm
- 9033511 2.5mm
- 9033503 3.0mm
- 9033510 3.5mm
- 9033504 4.0mm
- 9033505 5.0mm
- 9033506 6.0mm
- 9033508 8.0mm
- 9033525 Assorted
Five each:
2.0mm 3.0mm,
3.5mm, 4.0mm, 5.0mm



Box of 50

- 9033519 1.5mm
- 9033512 2.0mm
- 9033521 2.5mm
- 9033513 3.0mm
- 9033520 3.5mm
- 9033514 4.0mm
- 9033515 5.0mm
- 9033516 6.0mm
- 9033518 8.0mm

1.5mm



2mm



2.5mm



3mm



3.5mm



4mm



5mm



6mm



8mm



Dermabraders


Grit Selection



Regular



Coarse



Extra Coarse

"PBS" Diamond Dermabrasion (shown actual size)

All dermabraders fitted with $\frac{3}{32}$ " shaft size.



17 x 10

2026103 Regular

2026106 Coarse



17 x 8

2026112 Regular

2026115 Coarse

2026118 Extra Coarse



17 x 6

2026121 Regular

2026124 Coarse

2026127 Extra Coarse



17 x 4

2026130 Regular

2026133 Coarse

Cone



2026150 Regular

2026153 Coarse

2026156 Extra Coarse

Pear



2026159 Regular

2026162 Coarse

2026165 Extra Coarse

Cylinder



2026168 Regular

Forceps



Dressing • 9085276 (4 1/2")



Dressing • 9085278 (5 1/2")



Dressing • 9085277 (6")



Tissue • 9085282 (1x2, 4 1/2")



Tissue • 9085283 (1x2, 5 1/2")



Tissue • 9085284 (1x2, 6")



Iris Dressing - Straight • 9085300 (4")



Iris Dressing - Curved • 9085301 (4")



Iris Tissue - Straight • 9085303 (4")



Iris Tissue - Curved • 9085304 (4")

All in ages are actual size.

Instruments


Forceps



Adson Dressing • 9085295 (4 3/4")



Adson Tissue • 9085296 (1 x 2, 4 3/4")



Adson-Brown Tissue • 9085299 (4 3/4")

Hand Instruments



Probe-Packer • 1038020 (DE, 6")



Spatula-Packer • 1038021 (DE, 5 3/4")



Excavator-Packer • 1038022 (DE, 5 3/4")



Probe-Excavator • 1038023 (DE, 6")

Hand Instruments



Locke Elevator • 1038035 (Narrow)



Locke Elevator • 1038036 (Wide)



Freer Elevator • 1031260 (Sharp/Blunt, 7 1/2")