


Double-Ended Disposable Ear Curettes

Unmatched Performance and Versatility



Double-Ended *Increases Efficiency and Cost Savings*

Position Indicators™ *Facilitate Guidance and Positioning*


Disposable *Offers Safety and Convenience*

Flexible Tips *Allow for Comfortable Insertion*


Three Different Shapes *Meet the Needs of Any Practice*

Premier®'s **Double-Ended** Disposable Ear Curettes feature the functionality of metal curettes with the comfort and convenience of disposables. Available in three styles for all your procedural needs.


Mini-Loop #1038043
For hardened wax.



Curve-Loop #1038042
For general procedures.



Child-Scoop #1038041
For smaller children and infants.



Packaged in boxes of 50.
Contents are non-sterile.

Available from:

#8610620PPL50110M

