



Premier® Medical Products Company
1710 Romano Drive, P.O. Box 4500, Plymouth Meeting, PA 19462, USA
610-239-6000 • 888 PREMUSA (773-6871) • Fax: 1-610-239-6171
medinfo@premusa.com • www.premusa.com



Premier® Medical Products



Premier®
Medical Products



Dear Customers & Healthcare Professionals:

We at Premier are truly delighted in your interest in our products. Thank you for the time spent reviewing our new Medical Products catalog. This catalog, which details our entire line of medical products, represents a milestone in our growth. As a division of the family-owned Premier Dental Products Company, the catalog clearly announces Premier Medical's reputation as a provider of high quality, innovative products to the healthcare industry.

So even if you are familiar with Premier Medical, please take a few moments to browse the catalog. It is very likely you will be surprised with the various products offered to physician specialties other than those you have traditionally associated with Premier.

Your complete and total satisfaction is our only goal at Premier. We are very interested in hearing from you with regard to new products, new techniques or a critique of the catalog. All comments or suggestions are welcome on how we can better service you and more importantly, those you serve.

Sincerely,



Kevin J. Brown
Vice President
Premier Medical Products

*All Premier Medical Products Are Available Through
Your Medical or Surgical Supply Distributor*

Estimados Clientes y Profesionales:

Premier es encantado en su interés en nuestros productos. Gracias por tomar el tiempo para repasar nuestro nuevo catálogo de Productos Médicos. Este catálogo, que detalla nuestra línea completa de productos médicos, representa un suceso culminante en nuestro crecimiento. Como una división de Premier Dental Products Company, el catálogo claramente anuncia la reputación de Premier Medical como un proveedor de productos innovadores de alta calidad a la industria de salud.

Aunque usted sea familiar con Premier Medical, por favor tome unos momentos de hojear el catálogo. Es muy probable que usted será sorprendido por varios productos ofrecidos a especialidades mas otras que usted tradicionalmente se ha asociado con Preimer.

Su satisfacción completa y total es nuestro único objetivo en Premier. Estamos muy interesados en tener noticias de usted con respeto a productos nuevos, técnicas nuevas o una crítica del catálogo. Todos los comentarios o sugerencias son bienvenidos sobre como podemos atenderle mejor, y más importante, atender aquellos que usted sirve.

Sinceramente,



Kevin J. Brown
Vicepresidente
Premier Medical Products

*Todos los Productos Médicos de Premier son Disponibles
Por Su Distribuidor de Suministro Médico o Quirúrgico*



All U.S. and Canadian inquiries should be directed to:

Premier® Medical Products

1710 Romano Drive, P.O. Box 4500, Plymouth Meeting, PA 19462, USA
1-888-773-6872 (U.S. & Canada) • 1-610-239-6000 • Fax: 1-610-239-6171
E-mail: medinfo@premusa.com • www.premusa.com

All Latin American and Overseas inquiries should be directed to:

Premier® Dental Products Company

International Division

1710 Romano Drive, P.O. Box 4500, Plymouth Meeting, PA 19462, USA
610-239-6000 • Fax: 610-239-6171 • E-mail: intlinfo@premusa.com

2001e Electrosurgery Generator2, 27

A

Accu-Mold™47
Air Filtration System (Smoke Evacuator).....3, 27
Air-Lon™
Tracheostomy and Laryngectomy Tubes42

B

Biopsy Forceps (Gynecology)11
Biopsy Punches, Skin, Disposable (Uni-Punch) ...29
Bipolar Forceps (Electrosurgery)28
Brite Shield™ (Instrument Cleaner)59
Burs, Podiatry46

C

Cannula, Uterine (and replacement parts)22
Cart4
Cone Biopsy Electrodes6
Cryosurgery (Nitrospray Plus)25
Curettes, Endocervical19
Curettes, Endometrial18
Curettes, Uterine17

D

Decannulation Stoppers, Tracheostomy43
Dermabraders30
Dilators, Cervical16
Dilators, Uterine14-15
Diamonds, Dermabraders30
Diamonds, Podiatry46

E

Ear Curettes (Double-Ended, Disposable)37
Electrodes, Cone Biopsy (TruCone™)6
Electrodes, Cutaneous27
Electrodes, LOOP, Disposable7
Electrodes, LOOP, Reusable7
Electrosurgery Generator, 2001e2, 27
Elevators (Locke, Freer)35

F

Foot Controller, Electrosurgery5
Forceps (Tissue and Dressing)31-32

G

Glass Bead Mirror Warmer38

H

Hand Instruments31-33
Hemostatic Agent (Lumicain®)47
Hemostatic Agent (Monsel's)9

I

Instrument Cleaner and Care59
Iodine Solution, Strong (Lugol's)9

L

Laryngectomy Tubes, Metal41
Laryngectomy Tubes, Nylon (Air-Lon™)42
Lateral Vaginal Retractor8, 13
"LOOP" Electrodes6-7
"LOOP" Surgery Accessories3-10
"LOOP" Surgery Instruments5, 8
Lugol's Solution9
Lumicain®47

M

Medicaments, "LOOP" Surgery9
Moldable Silicone Compound (Accu-Mold™)47
Monsel's Solution - Thickened9

N

Nail Instruments (Nipper)46-47
Needle Adapters10
Needle, Cervical Block, Reinforced10
Needle Holders52-54
Nitrospray Plus25
Nucleus Knife34

O

Os Dilators16

P

PBS Diamond Burs44
Pessary Products10-11

R

Retractors8, 13, 44, 57

S

Salinocaine47
Scissors, Surgical48, 50, 55-56
Sheet Holder, Bethea23
Skin Biopsy Punch, Disposable (Uni-Punch)29
Smoke Evacuator3, 27
Sounds, Uterine17
Specula, Endocervical13
Specula, Vaginal8, 12-13
Suction Tubes58
Surgical Headband59

T

Tenacula8, 20-21
Tongue Depressors44
Tracheostomy Tubes, Accessories43
Tracheostomy Tubes, Metal39-40
Tracheostomy Tubes, Nylon (Air-Lon™)42
TruCone™6

U

Uni-Punch®29

V

Verruca (Wart) Removal (Salinocaine)47



8 ML NDC 48783-112-08
premier Monsele's
(Ferric Sub sulfate Solution)
Caution: Restricted to use by or on
the order of a physician.

2001e
Electrosurgical Generator

120

80

Power Control

Return Electrode Alarm

Monopolar

Bipolar

Coagulation

Fulguration

Splicer

premier

Power

Return Electrode

Pool Control

Monopolar

Bipolar

Gynecology

[“LOOP” Surgery Products](#)

[Medicaments](#)

[Pessary Products](#)

[Biopsy Forceps](#)


[Specula](#)


[Uterine Dilators](#)

[Cervical Dilators](#)

[Uterine Sounds](#)

[Instruments](#)





*Designed with
the alternate site
in mind!*

Model 2001e Electrosurgical Unit

Full-Featured Electrocauter
for the Office

- 120 watts Maximum Power Output
Perform office-based procedures with confidence.
- Five Distinct Output Modes
Cut, Coagulation, Blend, Fulguration and Bipolar outputs provide unsurpassed versatility.
- Isolated Output Transformers and Patient Circuit Monitor
Ensure safety to both patient and staff.
- Flush Membrane Front Panel
No knobs or switches to adjust, easy to clean and disinfect.
- Digital Diagnostic Error Codes
Readily identify and protect against potential fault conditions.
- Hand Control with Optional Foot Control
Fully adaptable to almost any procedure or electrode.
- Meets the Following Safety and Performance Standards:
CSA C22.2 NO 601.1-M90, UL 2601-1-UL, IEC 60601-2-2 (1998-90)

9006501 Model 2001e Electrocauter Unit (110/240VAC)
(Includes: hospital-grade power cords (110V and 220V) disposable hand-control pencil, disposable patient pad, clamp cord assembly and manual)

SPECIFICATIONS:

Maximum Power Outputs by Operating Mode:
 Pure Cut = 120 Watts @ 500Ω, 357kHz ± 50 kHz
 Blend = 90 Watts @ 800Ω, 357kHz ± 50 kHz
 Coagulation = 80 watts @ 1000Ω, 475kHz ± 19 kHz
 Fulguration = 40 watts @ 1000Ω, 410kHz ± 50 kHz
 Bipolar = 30 watts @ 200Ω, 535kHz ± 29 kHz

AC Input Current: 110/240VAC
 Mains line frequency range (nominal): 50-60Hz
 Power consumption: 540VA

Dimensions:
 Size – 10 1/4”(W) x 12”(L) x 6” (H)
 (26cm x 30.5cm x 15.2cm)
 Weight – 14 lbs (6.5kg)

Two (2) year warranty



Evacuator/Air Filtration System

The Premier Evacuator/Air Filtration System with its distinctly quiet operation (60 decibels at 1.0 meter) is suitable for electrosurgical or laser procedures. ‘Responsible’ filter pack conveniently and safely collects plume (smoke and particulate matter to 0.2 microns at 99.999% efficacy¹) and neutralizes odors. Each filter pack comes with 3/8” tubing to attach to specula fitted with a smoke-removal tube.

1. Filter rating is for 15 procedures or 45 minutes of accumulated filtration time, whichever comes first.

9006300 Premier Evacuator (110V)

Includes One (1) Filter

- 9006298 Premier Evacuator (220V)
- 9006290 Replacement Filter Packs (5 ea/case)
- 9006292 Reusable Smoke Tube (for accessible sites)
- 9006296 Evacuator Footswitch (not shown)

SPECIFICATIONS:

AC Input Current:

- 9006300 - 110/120VAC, 60Hz, 15 amps
- 9006298 - 220/240VAC, 50Hz, 8 amps

Dimensions:

- Size - 10”(W) x 10”(L) x 14”(H)
(25.4cm x 25.4cm x 35.6cm)
- Weight - 20lbs. (9.1kg)

Two (2) year warranty

“LOOP” Surgery Products

ESU/Evacuator Procedure Cart

- Constructed of durable, non-conductive ABS plastic
- 300 lb. (135kg) capacity
- Three large, spacious shelves
- Transport handle and quiet, directional casters included
- Snap-together assembly (no tools required)

9006560 Premier ESU/Evacuator Procedure Cart

Dimensions:

Shelf Size - 16" x 24"
 (40.64cm x 60.96cm)
 Height - 34 7/8" (88.6cm)
 Weight - 32 lbs. (14.4kg)
 Capacity - 300 lbs. (135kg)



Optional procedure cart pulls the whole package together!

Hand-Control Pencil



Activates generator's cut and coagulation output modes. Fitted with a 10' cord, accepts standard shaft size (3/32") electrodes. (Compatible with other manufacturers' generators.)

9006505 Hand-Control Pencil, Disposable, 5/pack

9006515 Hand-Control Pencil, Reusable, one each

Foot-Control Pencil



For use with Foot Controller (9006525). Ten foot cord comes with standard banana plug (Bovie plug adapter needed for use with 2001e). Accepts standard shaft size (3/32") electrodes. (Also can be used with other manufacturers' generators.)

9006530 Foot-Control Pencil, Disposable, 5/pack

9006585 Bovie Plug Adapter, one each

Dual-Foot Controller



Dual operating mode conveniently selects either cut or coagulation operating modes. Standard color-coded pedals and fluid resistant cover (not shown) included. Also for use with bipolar accessories (see page 28).

9006525 Dual-Foot Controller

(For use with the Premier 2001e generator only.)

"LOOP" Surgery Products

Disposable Dispersive Patient Pads



9006540 Dispersive Pad², Single Cell, Disposable, 5/pack



9006550 Dispersive Pad³, Dual Cell, Disposable, 5/pack
For use with pad monitoring systems.

9006535 Replacement Clamp Cord Assembly (not shown)
For connecting unrecorded patient plates to the Premier 2001e generator.

PolyPlate™

Reusable Dispersive Patient Plate



Innovative - Utilizes advanced material technology.

Safe - Eliminates the risk of skin rashes, allergic reactions and possible lesions that can occur with the use of adhesive pads or pads requiring gels.

Economical - Saves costs, easy to clean and reduces medical waste.

General usage guidelines are printed on the non-patient contact side of the plate. Package insert provides detailed product instructions.

For use with all electrosurgery generators without a dispersive plate monitoring system. Compatible with Premier and non-Premier generators.

9006545 Reusable PolyPlate™, Single Cell

2. All Premier dispersive pads are unrecorded.
3. Not for use with 2001e.

“LOOP” Surgery Products



TruCone™



Rotational Cone Biopsy Electrode

- Easy to Use
No awkward hand movements required.
- Physician Assured
Stabilizing prongs increase procedure confidence.
- Improved Diagnosis
Produces a clean, smooth cut with minimum thermal artifact. Sample conveniently secured.
- Patient Friendly, Safe
Insulated components ensure safety to the patient, physician and staff.
- Eight (8) Different Sizes
The largest number of sizes available.

Packaged five per box. Each TruCone is packaged with instructional pathological guidelines.

	Diameter	Depth
9006110	18mm	12mm
9006115	18mm	15mm
9006125	24mm	12mm
9006130	24mm	15mm
9006140	30mm	12mm
9006145	30mm	15mm
9006165	30mm	9mm
9006170*	12mm	21mm

9006175 Assorted Pack (3)
1 each: 18mm x 12mm, 18mm x 15mm, 24mm x 15mm



Obtain Precise Cervical Cone Biopsies – Effortlessly! *It is as easy as...*

1

2

3

With the electro-surgical generator in the active mode, insert the appropriately sized TruCone™ into the endocervical canal. (Use of a foot-controlled generator recommended.)

After the stabilizing prongs have secured the cervix, gently rotate the shaft spindle for one complete rotation. Hub serves as a rotation guide.

Deactivate generator and remove TruCone instrument containing secured sample from vaginal cavity.

"LOOP" Electrodes

All Premier Electrodes are made with tungsten wire and fitted with a standard $\frac{3}{32}$ " diameter, 5" long shaft. Compatible with most electrosurgical generators on the market.


Reusable "LOOP" Electrodes

5 to 6 uses with proper care. Available individually.



Disposable "LOOP" Electrodes

Sterile, single-use electrodes. Available in 5-packs.



“LOOP” Surgery Products

“LOOP” Surgery Instruments

Designed specifically for use during LOOP and other electrosurgical procedures⁴, the non-conductive coating on these instruments ensures the safety of both patient and physician.

- Scratch Resistant, Non-conductive Coating
- Fully Autoclavable

Use and Care Instructions

- Can be autoclaved or cold soaked
- Do not use with ultrasonic cleaner
- Do not use if coating is compromised
- Do not use with dry heat sterilization

Coated Graves Specula

All Premier “Loop” Surgery specula are affixed with a smoke-removal tube. The tube, mounted on the upper blade, provides safe and convenient evacuation⁵ of electrosurgical plume from the vaginal cavity.

9086165 Wide Medium (1 1/4" x 4 1/2" blades)

9086166 Wide Large (1 3/8" x 5" blades)

9086500 Small (3/4" x 3" blades)

9086501 Medium (1 1/4" x 4 1/2" blades)

9086502 Large (1 3/8" x 5" blades)



Wide Large
(not actual size)

Coated Lateral Vaginal Retractors



With Offset Blades
(not actual size)

9086186 Lateral Retractor
with offset blades

Fully opened, this retractor provides an additional 1" of retraction.

Overall Length - 10"
Blade Size - 1" x 2 1/2"
Fully opened dilation - 3 1/2"


9086182 Lateral Retractor
(not shown)

Overall Length - 10"
Blade Size - 1" x 2 1/2"
Fully opened dilation - 2 1/2"

Coated Hooks (not actual size)

9086515
Iris Hook -
10"

9086520
Tenaculum
Hook - 9 1/2"



Coated Tenaculum (not actual size)



9086525 Schroeder
Tenaculum - 9 3/4"


4. These instruments are not recommended for use with surgical lasers.
5. Requires use of an air filtration system.

Medicaments

Monseľ's Solution



This thickened formulation of the popular hemostatic agent can be applied directly to the surgical site following excision of diseased tissue.



9045055 Box of twelve (12) 8ml bottles with applicators

Lugol's Solution



Lugol's is a strong iodine contrast medium for identification of diseased cervical tissue.



9045056 Box of twelve (12) 8ml bottles with applicators

Popular in LOOP, Biopsy and Colposcopic Procedures

Thickened formulation "adheres" to the surgical site for effective control of post-op bleeding.




Pessary Products

Pessary Products

For non-surgical management of uterine prolapse, cystocele and rectocele

- 100% silicone
- Easy to insert, remove and clean
- Single-patient use

Ring Pessary - This very common pessary most often is used for first- to second-degree prolapse. Ring Pessaries, available in nine sizes, are configured without and with support. Ring Pessaries with support also can be used on an accompanying cystocele. Application is accomplished easily simply by folding the ring before insertion.



Size	Diameter (inches)	Ring w/o Support Item #	Ring w/ Support Item #
#1	2.00	1040111	1040101
#2	2.25	1040112	1040102
#3	2.50	1040113	1040103
#4	2.75	1040114	1040104
#5	3.00	1040115	1040105
#6	3.25	1040116	1040106
#7	3.50	1040117	1040107
#8	3.75	1040118	1040108
#9	4.00	1040119	1040109

Cube Pessary - Designed for second- or third-degree prolapse or procidentia, as well as a cystocele and rectocele. Cube Pessaries feature drainage holes. Additionally, a silicone tie affixed to each cube aids in removal. Cube Pessaries are available in 10 sizes. Sizes range from 2.5cm (1") O.D. to 7.5cm (3") O.D.



Size	Width (cm)	Cube w/ Holes Item #
#0	2.5	1040400
#1	2.9	1040401
#2	3.2	1040402
#3	3.7	1040403
#4	4.1	1040404
#5	4.5	1040405
#6	5.5	1040406
#7	6.5	1040407
#8	7.0	1040408
#9	7.5	1040409

Gellhorn Pessary - Used for a second- to third-degree prolapse or procidentia. With drainage holes in the base, Gellhorn Pessaries are available in nine sizes. The knob of the gellhorn easily folds over for insertion and when properly inserted, rests on the posterior vaginal wall to provide support.



Size	Diameter (inches)	Gellhorn Item #
#0	1.50	1040200
#1	1.75	1040201
#2	2.00	1040202
#3	2.25	1040203
#4	2.50	1040204
#5	2.75	1040205
#6	3.00	1040206
#7	3.25	1040207
#8	3.50	1040208

Donut Pessary - For use with third-degree prolapse or procidentia as well as cystocele and rectocele. Donut Pessaries, while having a substantial amount of material between the inside and outside diameter, are made of soft material that is compressed easily to aid in insertion. The donut is available in seven different sizes.



Size	Diameter (inches)	Donut Item #
#0	2.00	1040300
#1	2.25	1040301
#2	2.50	1040302
#3	2.75	1040303
#4	3.00	1040304
#5	3.25	1040305
#6	3.50	1040306



1040600 **Pessary Fitting Set** - Set of six different rings (#2 through #7) used in pessary size selection for individual patients. The Fitting Set comes with a cross-reference chart that translates ring size into corresponding size of other pessary styles.

Needle Adapters

Reusable Stainless Steel Needle Adapters

These needle adapters allow for economical extension of an ordinary hypodermic syringe. The now “extended” needle easily reaches the paracervical and uterosacral areas to provide improved control for local infiltration of anesthetic agents. Slightly malleable for additional convenience. Four different lengths. Luer lock connectors.



3" Length
• 1032453



4" Length
• 1032454



5" Length
• 1032455



6" Length
• 1032456


Biopsy Forceps

German Biopsy Punch Forceps

German Precision Instruments

- Exact Bite for Clean Tissue Extraction
- Firm Hold of Biopsy Sample for Confident Procedures
- Spring Grip for Optimal Handling
- Swivel Lock to Ensure Sample Retention
- Popular 10" Length
- Outstanding Value, One-Year Unconditional Warranty

Made in Germany from high-grade surgical stainless steel, Premier®'s Biopsy Punch Forceps offer optimal performance and ultimate quality. Through superior craftsmanship, these four classic shapes afford the highest level of clinical care.



9086420 Tischler

This clinical standard features a unique bite mechanism which grasps tissue and holds it securely during excision.

Bite: 7mm x 3mm



9086430 Baby Tischler

A mini or Baby Tischler provides a smaller sample. Its rounded bite minimizes patient discomfort and reduces the occurrence of bleeding.

Bite: 4.2mm x 2.3mm



9086410 Burke

This small-toothed forceps allows for sample retrieval with minimal tissue penetration.

Bite: 5mm x 3mm



9086440 Kevorkian

This classic pattern offers the largest available jaw and the serrated teeth afford strong grasping ability.

Bite: 9mm x 3mm

Specula



**Graves Specula
Upper Blade**



**Pederson Specula
Upper Blade**




Graves

9086100 Small (3/4" x 3" blades)

9086101 Medium (1 1/4" x 4 1/2" blades)

Medium
(not actual size)



Pederson


9086110 Small (3/4" x 3" blades)

9086111 Medium (1 1/4" x 4 1/2" blades)

9086112 Large (1 3/8" x 5" blades)

Medium
(not actual size)

Specula



Weisman-Graves

9086140 LLO (1 1/4" x 4 1/2" blades)

9086141 RLO (1 1/4" x 4 1/2" blades) (not actual size)



Weisman-Graves
Specula
Upper Blade


Diameters, French Scale

	mm	inch
○ 9	3.0	.118
○ 10	3.33	.131
○ 11	3.67	.144
○ 12	4.0	.158
○ 13	4.33	.170
○ 14	4.67	.184
○ 15	5.0	.197
○ 16	5.33	.210
○ 17	5.67	.223
○ 18	6.0	.236
○ 19	6.33	.249
○ 20	6.67	.263
○ 21	7.0	.275
○ 23	7.67	.302
○ 25	8.33	.328
○ 27	9.0	.354
○ 29	9.67	.380
○ 31	10.33	.407
○ 33	11.0	.433

Uterine Dilators

Hank French scale (not actual size)

Set of Six • 1030541



9/10 • 1030530

11/12 • 1030532

13/14 • 1030534

15/16 • 1030536

17/18 • 1030538

19/20 • 1030540

Pratt French scale (not actual size)

Set of Eight • 1030580



13/15 • 1030565

17/19 • 1030567

21/23 • 1030569

25/27 • 1030571

29/31 • 1030573

33/35 • 1030575

37/39 • 1030577

41/43 • 1030579

French	35	37	39	41	43
mm	11.67	12.33	13.0	13.67	14.33
inch	.459	.485	.512	.538	.564

Uterine Dilators

Hegar (not actual size) **Set of Eight** • 9030559



3/4 (3mm/4mm) • 9030544



5/6 (5mm/6mm) • 9030546



7/8 (7mm/8mm) • 9030548



9/10 (9mm/10mm) • 9030550



11/12 (11mm/12mm) • 9030552




13/14 (13mm/14mm) • 9030554



15/16 (15mm/16mm) • 9030556



17/18 (17mm/18mm) • 9030558




Kleegman (not actual size)

Stainless 1.8mm/2.5mm • 1030581



PTFE 1.5mm/2.1mm • 1030582

Diameters, Millimeters

○ 1  11

○ 2  12

○ 3  13

○ 4  14

○ 5  15

○ 6  16

○ 7  17

○ 8  18

○ 9

○ 10

Cervical Dilators

Each Premier® Os dilator is manufactured to exacting standards. Made of PTFE (Polytetrafluoroethylene), which is malleable, 'warm' and self-lubricating, these instruments deliver excellent patient comfort. Usage and care instructions are provided below.

Os Locator is distinguished with a single ring and a yellow band around the handle. It is particularly useful in locating the external Os in patients presenting cervical stenosis.


Canal Finder has a double ring and a blue band around the handle. The slightly longer tip insinuates the cervical canal, while producing dilation.

Fundus Dilator is fitted with a triple ring and a green band on its handle. Having the longest tip of the set, this particular instrument can be passed through the entire length of the cervical canal and will open the internal Os.

The following guidelines ensure proper care and performance of the Os Dilators:

- Do not bend the instrument more than 90° from its longitudinal axis.
- Inspect instrument for cracks or breakage prior to EACH use.
- Clean and sterilize before each use.
- Os Dilators can be cleaned and sterilized by any approved method. (Dry heat sterilization not recommended.)

OS Dilators Reusable (not actual size) **Set of Three** • 1030585



OS Locator
• 1030586
(Yellow Ring)

Canal Finder
• 1030587
(Blue Ring)

Fundus Dilator
• 1030588
(Green Ring)

Uterine Sounds

Sims



Graduated/cm • 1030522



Graduated/cm (2mm) • 1030524



Sim's Handle Style (not actual size)
 All tips shown actual size.
 All Premier Sounds are silver plated and malleable.

Endometrial Biopsy Curettes

2mm (not actual size) All tips shown 200% actual size.


Novak



2mm • 1030312

3mm • 1030313

4mm • 1030314



Endocervical Curettes

Kevorkian-Young

(not actual size)



Malleable • 1030330
(not shown)

Malleable with Basket • 1030331
(not shown)



Rigid • 1030332



Rigid with Basket • 1030333

Duncan


(not actual size)




Rigid • 1030343

Instruments


Atraumatic Jaws



Tenaculum (not actual size)




Cervical Tenaculum • 9030610
(1/4" with Atraumatic Jaws)




Cervical Tenaculum • 9030611
(1/8" with Atraumatic Jaws)



Emmett Tenaculum Hook #1
• 1030603 (9", Half Curve)




Schroeder Uterine Tenaculum
 • 9030613 (10")



Hulka Uterine Tenaculum
 • 1030740 (11 1/4" with Sound)



Hulka Uterine Tenaculum
 • 1030742 (11 1/4" with Insulated Sound)



Hulka Uterine Tenaculum
 • 1030744 (11 1/4" with Atraumatic Jaw and Sound)

Instruments

Cannula Instruments

Jarcho Self-Retaining Cannula • 1030840
(not actual size)



Kahn Uterine Trigger Cannula
• 1030842 (with 3 Spare Tips)
(not actual size)



These instruments are particularly useful in the administration of radiographic dyes and other solutions into the uterine cavity. Flexible malleable tips prevent the need for sounding or dilation.

Replacement Parts for Cannula Instruments

1032100 Rubber Acorn (3/pack)

1032219 PTFE Tip, Short (3/pack)

1032220 PTFE Tip (3/pack)


1032221 PTFE Acorn

1032222 Collar/Screw (with Rubber Gasket)


1032223 PTFE Acorn for Corson

1032224 Collar Screw for Corson

Miscellaneous Instruments



Polyp Forceps, Narrow • 9030625 (10")
(not actual size)



IUD Removal Hook • 1030655
(not actual size)




IUD Removal Hook, Universal • 1030656
(not actual size)

Sheet Holders - OB/GYN

Bethea Sheet Holder
• 1050222 (1 Pair)
(not actual size)





Dermatology

[Liquid Nitrogen Cryosurgery Instruments](#)

[Tip Attachments/Miscellaneous](#)

[Electrosurgical Unit](#)

[Smoke Evacuator](#)

[Bipolar Accessories](#)

[Disposable Biopsy Punch](#)

[Dermabraders](#)

[Instruments](#)

Cryosurgery

Nitrospray Plus

Trigger Operated, Closed System,
Liquid Nitrogen Cryosurgical Unit

- Improved Seal⁶
Eliminates o-ring positioning.
- Grooved Head
Easy to use, comfortable to hold.
- Curved Trigger
Unobstructed view of treatment site.
- Notched Head Grip
Easy to remove/install cap.
- Three (3) year warranty

**30-Day
Unconditional
Warranty**

Static holding time:
8 - 12 hours

Intermittent use:
3 - 6 hours of patient treatment time



1006060 Nitrospray Plus
(16oz fill capacity)
Includes 5 spray tips,
tip protector and
base holder



1006065 Nitrospray Plus Lite
(10oz fill capacity)
Includes 5 spray tips,
tip protector and
base holder



6. Pressure seal ring retrofit kit (1006066) to eliminate o-ring on units prior to RP-11061 and LP-10952 is available. Not to be used with RP-10454 to RP-10561. Contact Customer Service for assistance.

Cryosurgery - Tip Attachments

Spray Tips

(shown actual size)



Spray Tip, 16 gauge • 1006516
(purple)



Spray Tip, 17 gauge • 1006517
(white)



Spray Tip, 18 gauge • 1006518
(pink)



Spray Tip, 19 gauge • 1006519
(brown)



Spray Tip, 20 gauge • 1006520
(yellow)

Spray Tips, Set of Five • 1006525
one of each size (16, 17, 18, 19, 20 gauge)

Cutaneous Probes*

(not actual size)



Cutaneous Probe, 2mm • 1006602



Cutaneous Probe, 4mm • 1006604



Cutaneous Probe, 6mm • 1006606



Cutaneous Probe, 8mm • 1006608



Cutaneous Probe, 10mm • 1006610



Cutaneous Probe, Sharp • 1006620

Cutaneous Probes, Set of Five • 1006615
(2mm, 4mm, 6mm, 8mm, 10mm)

*Reusable vent tube provided with each probe.

Cryocones

(not actual size)



Cryocone, 2.5mm • 1006625



Cryocone, 5mm • 1006635



Cryocone, 10mm • 1006640



Cryocone, 20mm • 1006650




Cryocone, 25mm • 1006665

Cryocones, Set of Five • 1006670
(2.5mm, 5mm, 10mm, 20mm, 25mm)

Miscellaneous

(not actual size)



Bent Spray Extension • 1006530



Spray Tip Protector • 1006527



Needle Extender, Rigid • 1006700



Brymill Accessory Adapter • 1006521

1006730 Replacement Bulb/Hose for Original Nitrospray

1006735 2" Needle Extender



**Model 2001e
Electrosurgical Unit**

Full-Featured Electrosurgery
for the Office

Designed with the alternate
site in mind!

Increased capabilities beyond
coagulation-only units!

9006501 Model 2001e Electrocauter Unit
(110/240VAC)
(Includes: hospital-grade
power cords (110V and 220V),
disposable hand-control
pencil, disposable patient pad,
clamp cord assembly and manual)

*See page 2 for a full description, specifications
and warranty information.*



**Evacuator/Air Filtration
System**

- 9006300 Premier Evacuator (110V)
Includes One (1) Filter

- 9006298 Premier Evacuator (220V)

- 9006290 Replacement Filter Packs (5 ea/case)

- 9006292 Reusable Smoke Tube
(for accessible sites)
Length: 8 7/8", Diameter: 3/4"

- 9006296 Evacuator Footswitch (not shown)

- 9006297 Hands-Free Tube Holder (not shown)
Eliminates the need of an assistant,
while providing a full range of motion.

*See page 3 for a full description, specifications and
warranty information.*

Electrosurgery

Bipolar Coagulating Forceps

Further expand the capability of the 2001e Electrosurgical Unit with Premier Bipolar Forceps. Fully autoclavable, with a standard two-prong connector, these forceps are compatible with other generators and bipolar cord sets.

Must be used with foot control (refer to 9006525 Dual Foot Controller below).



Because the output signal is contained between the forceps tips, no grounding pad is required to complete the electrical circuit.

This makes bipolar the preferred mode of coagulation when doing simple procedures around the head, face and upper chest.

As with any electrosurgical procedure, extreme caution should be used on patients with cardiac pacemakers, diagnostic monitoring systems or infusion pumps.

(Illustration above for descriptive purposes only. Not representative of actual procedure.)



9006308 5 1/2" Straight



9006310 7 1/2" Bayonet, Fine Tip



9006309 5" Curved



9006314 7" Straight, Fine Tip

Contact Premier Customer Service for the availability of additional forcep patterns.

9006313 Standard Bipolar Adapter Cord, reusable (Not Shown)
Banana plug male connectors are compatible with most manufacturers' generators. Accepts forceps with male 2 pin connectors.

Dual Foot Controller

Fluid resistant, dual operating mode conveniently selects either cut or coagulation operating modes. Standard color-coded pedals included.



9006525 Dual Foot Controller

(For use with the Premier 2001e generator only.)

Disposable Biopsy Punch



Seamless razor-sharp stainless steel blade for exact excisions.


Ribbed handle offers superior control.

Diameter reading for quick identification.

Hollow handle provides easy specimen retrieval.

Uni-Punch®

Sterile, Seamless, Razor Sharp Blades



Box of 25

- 9033509 1.5mm
- 9033502 2.0mm
- 9033511 2.5mm
- 9033503 3.0mm
- 9033510 3.5mm
- 9033504 4.0mm
- 9033505 5.0mm
- 9033506 6.0mm
- 9033508 8.0mm
- 9033525 Assorted
Five each:
2.0mm 3.0mm,
3.5mm, 4.0mm, 5.0mm



Box of 50

- 9033519 1.5mm
- 9033512 2.0mm
- 9033521 2.5mm
- 9033513 3.0mm
- 9033520 3.5mm
- 9033514 4.0mm
- 9033515 5.0mm
- 9033516 6.0mm
- 9033518 8.0mm

1.5mm



2mm



2.5mm



3mm



3.5mm



4mm



5mm



6mm




8mm



Dermabraders

Grit Selection



Regular



Coarse



Extra Coarse

"PBS" Diamond Dermabrasion *(shown actual size)*

All dermabraders fitted with $\frac{3}{32}$ " shaft size.



Forceps



Dressing • 9085276 (4 1/2")



Dressing • 9085278 (5 1/2")



Dressing • 9085277 (6")



Tissue • 9085282 (1x2, 4 1/2")



Tissue • 9085283 (1x2, 5 1/2")



Tissue • 9085284 (1x2, 6")



Iris Dressing - Straight • 9085300 (4")



Iris Dressing - Curved • 9085301 (4")



Iris Tissue - Straight • 9085303 (4")



Iris Tissue - Curved • 9085304 (4")

All images are actual size.

Instruments

Forceps



Adson Dressing • 9085295 (4 3/4")



Adson Tissue • 9085296 (1 x 2, 4 3/4")



Adson-Brown Tissue • 9085299 (4 3/4")

Hand Instruments



Probe-Packer • 1038020 (DE, 6")



Spatula-Packer • 1038021 (DE, 5 3/4")



Excavator-Packer • 1038022 (DE, 5 3/4")



Probe-Excavator • 1038023 (DE, 6")

Hand Instruments




Locke Elevator • 1038035 (Narrow)



Locke Elevator • 1038036 (Wide)



Freer Elevator • 1031260 (Sharp/Blunt, 7 1/2")



ENT

[Double-Ended Disposable Ear Cures](#)

[Glass Bead Mirror Warmer](#)

[Tracheostomy & Laryngectomy Tubes](#)

[Tongue Depressors](#)

Double-Ended Disposable Ear Curettes

Premier®'s Double-Ended Disposable Ear Curettes feature the functionality of metal curettes with the comfort, safety and convenience of disposables. Available in three styles for all your procedural needs.

- Double-Ended**
Increases Efficiency and Cost Savings
- Position Indicators™**
Facilitate Guidance and Positioning
- Disposable**
Offers Safety and Convenience
- Flexible Tips**
Allow for Comfortable Insertion
- Three Different Shapes**
Meet the Needs of Any Practice



Mini-Loop

Curve-Loop

Child-Scoop

Premier's unique Position Indicators™ provide accurate tip orientation in the ear canal.



Packaged in boxes of 50.
Contents are non-sterile.

- 1038043 Mini-Loop
(For hardened wax)
- 1038042 Curve-Loop
(For general procedures)
- 1038041 Child-Scoop
(For smaller children and infants)

Glass Bead Mirror Warmer



Unsurpassed convenience in warming mirrors⁷ during laryngeal examinations.

Warms mirrors to prevent fogging.

Meets the needs of the ENT clinician for improved laryngeal examinations... and offers these advantages!

- Hinged lid for convenient use.
- Rapid warm-up⁸ for greater efficiency.
- Thermal wall remains cool and safe to touch.
- Thermostatically controlled to maintain required temperature.
- UL approved.
- Dimensions: 3 1/4" x 5 5/8" x 3 1/4"

Simply immerse the mirror into the hot beads for three to five seconds when needed.

Cautions:

- Allow mirror sufficient time for cooling before use.
- Soiled or contaminated beads must be changed.

See product information sheet enclosed with unit for complete details.

1006131 Mirror Warmer and Glass Beads (110v)

1006132 Mirror Warmer and Glass Beads (220v)

1006105 Replacement Beads

Three (3) year warranty

7. Not approved for use as a sterilizer.

8. Inside well reaches a temperature of 250°F after 20 minutes.

Tracheostomy Tubes

Jackson Original

- Have a 90° curve and a swivel lock to hold the inner cannula in place.
- Neckplate Shapes
 Sizes 00-8 have a small, rectangular shape while sizes 9-10 have the larger oval-shaped neckplate.



Original size 6 shown

Regular Length

Item #	Size	MM OD	MM ID	MM Length	Ins.
1034300	00	3.4	1.6	33	1.30
1034310	0	4.6	2.8	38	1.50
1034311	1	5.2	3.5	44	1.73
1034312	2	6.1	4.0	51	2.01
1034313	3	7.1	5.1	58	2.28
1034314	4	8.2	5.8	62	2.44
1034315	5	9.2	6.8	68	2.68
1034316	6	10.2	7.7	69	2.72
1034317	7	11.2	8.7	69	2.72
1034318	8	12.2	9.6	69	2.72
1034319	9	13.3	10.2	69	2.72
1034320	10	14.3	11.1	69	2.72

Short Length

Follows the original Jackson curve and meets the need for lengths between the regular Jackson and the short Luer Tube.

Item #	Size	MM OD	MM ID	MM Length	Ins.
1034600	00	3.4	1.6	33	1.30
1034610	0	4.6	2.8	36	1.42
1034611	1	5.2	3.5	40	1.57
1034612	2	6.1	4.0	44	1.73
1034613	3	7.1	5.1	48	1.89
1034614	4	8.2	5.8	52	2.05
1034615	5	9.2	6.8	56	2.20
1034616	6	10.2	7.7	56	2.20
1034617	7	11.2	8.6	56	2.20
1034618	8	12.2	9.6	56	2.20
1034619	9	13.3	10.2	56	2.20
1034620	10	14.3	11.1	56	2.20

Construction

The complete line of Premier Tracheostomy and Laryngectomy tubes are manufactured of stainless steel to exacting standards and under strict quality controls.

Precise tolerances and high quality guarantee a number of important benefits:

Interchangeable

An interchangeable inner tube correctly fitted to the outer tube provides maximum airway and more efficient processing, cleaning and assembly. All Jackson Original and Improved Tracheostomy Tubes feature interchangeable construction of the same size.

Pilot

A pilot (obturator) is shaped to be flush with the outer cannula to facilitate insertion.

Dimensions

Dimensions include tube size, outside diameter, inside diameter and length. Measurements are taken along the inside curve for the tube length.

Tracheostomy Tubes



Improved size 7 shown

Jackson Improved

Regular Length

- Have a 90° curve and a rotating lock to hold the inner cannula in place.
- Neckplate Shapes
 Sizes 00-6 have a small, rectangular shape while sizes 7-10 have the larger oval-shaped neckplate.

Item #	Size	MM OD	MM ID	MM Length	Ins.
1034900	00	3.4	1.6	33	1.30
1034910	0	4.6	2.8	38	1.50
1034911	1	5.2	3.5	44	1.73
1034912	2	6.1	4.0	51	2.01
1034913	3	7.1	5.1	58	2.28
1034914	4	8.2	5.8	62	2.44
1034915	5	9.2	6.8	68	2.68
1034916	6	10.2	7.7	69	2.72
1034917	7	11.2	8.7	69	2.72
1034918	8	12.2	9.6	69	2.72
1034919	9	13.3	10.2	69	2.72
1034920	10	14.3	11.1	69	2.72

Short Length

Follows the regular Jackson curve and meets the need for lengths between the regular Jackson and the short Luer tube.

Item #	Size	MM OD	MM ID	MM Length	Ins.
1034930	00	3.4	1.6	33	1.30
1034940	0	4.6	2.8	36	1.42
1034941	1	5.2	3.5	40	1.57
1034942	2	6.1	4.0	44	1.73
1034943	3	7.1	5.1	48	1.89
1034944	4	8.2	5.8	52	2.05
1034945	5	9.2	6.8	56	2.20
1034946	6	10.2	7.7	56	2.20
1034947	7	11.2	8.7	56	2.20
1034948	8	12.2	9.6	56	2.20
1034949	9	13.3	10.2	56	2.20
1034950	10	14.3	11.1	56	2.20

Extra Long

Are 16mm (5/8 in.) longer than the standard Jackson in sizes 5 through 10. In order to maintain the correct arc they are not as long in size 4.

Item #	Size	MM OD	MM ID	MM Length	Ins.
1034964	4	8.2	5.8	72	2.83
1034965	5	9.2	6.8	85	3.35
1034966	6	10.2	7.7	85	3.35
1034967	7	11.2	8.7	85	3.35
1034968	8	12.2	9.6	85	3.35
1034969	9	13.3	10.2	85	3.35
1034970	10	14.3	11.1	85	3.35

Jackson Improved with Adapter

Regular Length

Are modified with a permanent 15mm adapter for use with positive pressure and anesthesia units.

Item #	Size	MM OD	MM ID	MM Length	Ins.
1034924	4	8.2	5.8	62	2.44
1034925	5	9.2	6.8	68	2.68
1034926	6	10.2	7.7	69	2.72
1034927	7	11.2	8.7	69	2.72
1034928	8	12.3	9.6	69	2.72



Improved size 8 shown

(Short or extra-long tubes may be fitted with a permanent adaptor.)

Laryngectomy Tubes

Jackson Original

- These tubes have a standard curve with a swivel lock on the neckplate.
- Neckplate Shapes
 Sizes 4-8 have a small, rectangular neckplate while sizes 9-10 have the larger oval-shaped neckplate.

Item #	Size	MM		Length	Ins.
		OD	ID		
1037214	4	8.2	5.8	50	1.97
1037215	5	9.2	6.8	53	2.09
1037216	6	10.2	7.7	55	2.17
1037217	7	11.2	8.7	55	2.17
1037218	8	12.2	9.6	55	2.17
1037219	9	13.3	10.2	55	2.17
1037220	10	14.3	11.1	55	2.17



Original size 6 shown

Jackson Improved

- These tubes have a standard curve with a rotating lock on the neckplate.
- Neckplate Shapes
 Sizes 4-6 have a small, rectangular neckplate while sizes 7-10 have the larger oval-shaped neckplate.

Item #	Size	MM		Length	Ins.
		OD	ID		
1037414	4	8.2	5.8	50	1.97
1037415	5	9.2	6.8	53	2.09
1037416	6	10.2	7.7	55	2.17
1037417	7	11.2	8.7	55	2.17
1037418	8	12.2	9.6	55	2.17
1037419	9	13.3	10.2	55	2.17
1037420	10	14.3	11.1	55	2.17



Improved size 7 shown

Martin

- The short length of this tube is standard on each of the available sizes and is attached to the plate at an acute downward angle.
- All Martin Tubes have the large, oval-shaped neckplate.

Item #	Size	MM		Length	Ins.
		OD	ID		
1036128	8	12.3	9.5	42	1.65
1036130	10	14.1	11.0	42	1.65



Air-Lon™ Tracheostomy and Laryngectomy Tubes

All Premier Air-Lon tubes include inner cannula, outer cannula and obturator. Length measured on outside curve from distal end to neck plate. Inner cannula extends about 1mm beyond tube end. Inhalation sets have the same dimensions as tubes and come with a 15mm adapter permanently attached to the inner cannula. Compatibility with Jackson-scale accessories (decannulation stoppers, etc.) not guaranteed.

- **Convenience** Lightweight, fully autoclavable tube made from an exclusive nylon material. Since Air-Lon tubes are non-metal, no 'warming' is required for patient comfort.
- **Durability** Air-Lon tubes resist scratching, denting and tarnishing. Friction-fitting components eliminate the need for locks or hinges that become damaged or break with use.
- **Sterilization Guidelines** May be boiled, autoclaved, gas sterilized or immersed in cold germicidal solutions (except those with phenol or formic acid base). Dry heat sterilization not recommended.

Tracheostomy Tubes



Tubes		O.D.	I.D.	Length	Length	Inhalation
Item #	Size	(mm)	(mm)	(mm)	(in.)	Sets
1050150	3	7	4	58	2.28	1050180
1050152	4	8	5	62	2.44	1050182
1050154	5	9	6	65	2.56	1050184
1050156	6	10	7	72	2.83	1050186
1050158	7	11	8	76	2.99	1050188
1050160	8	12	9	81	3.19	1050190

90° Curvature

Inhalation Sets with 15mm adapter



Laryngectomy Tubes



Tubes		O.D.	I.D.	Length	Length	Inhalation
Item #	Size	(mm)	(mm)	(mm)	(in.)	Sets
1050170	8	12	9	42	1.65	1050200
1050172	10	14	11	40	1.55	1050202

Martin Style Curve and Length

Inhalation Sets with 15mm adapter



Tracheostomy & Laryngectomy Accessories



Adapters

Jackson Improved Tracheostomy Tubes (regular length-adult sizes 4-8) are available as line items with a 15mm adapter. Short or extra long Improved tubes also may be fitted with a permanent adapter if required.



Fenestrations

Tubes may be fenestrated on one or both cannulae. Fenestrations are usually at a distance 1/3 the length of the tube-measured from the neckplate - and must be physician ordered. Fenestrated tubes are sometimes described as "speaking" tubes.

Penn Aspiring Tube 1050210 (not shown)

Tracheal Tube Brushes (not actual size)



Packaged in boxes of 12.
Contents are non-sterile.

9037702 Small (Tubes 1-6)

9037703 Large (Tubes 7-10)



Decannulation Stoppers

These color-coded plugs are available in four closure sizes in order to obstruct or partially obstruct the lumen of the tracheostomy tube during decannulation. Each stopper has a patient neck cord attached.

Yellow Stopper

(1/4 Closure)



Tube Size
0 - 1037800

1 - 1037801

2 - 1037802

3 - 1037803

4 - 1037804

Green Stopper

(1/2 Closure)



Tube Size
0 - 1037810

1 - 1037811

2 - 1037812

3 - 1037813

4 - 1037814

Purple Stopper

(3/4 Closure)



Tube Size
0 - 1037820

1 - 1037821

2 - 1037822

3 - 1037823

4 - 1037824

White Stopper

(Full Closure)



Tube Size
0 - 1037830

1 - 1037831

2 - 1037832

3 - 1037833

4 - 1037834

Tube Size
5 - 1037805

6 - 1037806

7 - 1037807

8 - 1037808

9 - 1037809

Tube Size
5 - 1037815

6 - 1037816

7 - 1037817

8 - 1037818

9 - 1037819

Tube Size
5 - 1037825

6 - 1037826

7 - 1037827

8 - 1037828

9 - 1037829

Tube Size
5 - 1037835

6 - 1037836

7 - 1037837

8 - 1037838

9 - 1037839

Tongue Depressors




Weder Retractor • 1050130 (Small)
(not actual size)



Weder Retractor • 1050132 (Large)
(not actual size)

Podiatry



["PBS" Diamond Burs](#)

[Accu-Mold™ - Moldable Silicone Compound](#)

[Medicaments](#)

[Nail Instruments](#)

"PBS" Diamond Burs - Nail Reduction & Tissue Debridement

Professional Treatment System For Hyperkeratotic Tissue And Hypertrophic Nails. Through the exclusive bonding process, which permanently bonds diamond to the bur surface, Premier burs offer unsurpassed product life and performance! More efficient than steel or carbide instruments and with greater patient acceptance because of reduced heat generation.

All diamond burs fitted with $\frac{3}{32}$ " shaft size.



Barrel • 2026005
(To reduce thick, deformed nails.)



Taper • 2026030
(For the treatment of incurvated nails or the nail groove.)



Pear • 2026010
(To smooth the rough edges of nails.)



Umbrella • 2026035
(For the reduction of hyperkeratoses, tylomas and heel callus.)




Bud • 2026015
(The most popular pattern for nail reduction.)



Small Olive • 2026040
(For working smaller areas.)



Short • 2026020
(For better control in treating smaller digits.)



Small Cylinder • 2026045
(For gross reduction of large nail areas.)



Olive • 2026025
(Similar to the PEAR and BUD.)



Diamond Bur System Kit • 2026050
A convenient kit of four nail and callus skin burs plus technique information. Includes one each:
Bud, Short, Taper, Umbrella



Clean-A-Diamond • 2014001
Dressing stone for diamond instruments.

Contact Premier Customer Service for bur cleaning and care guidelines.

Accu-Mold™ - Moldable Silicone Compound



9045050 Accu-Mold - 200 Gram Kit (yields approximately 15 devices)

Legendary in its convenience, Accu-Mold™ remains a significant advance over other moldable compounds. Simple to mix, easy to mold, Accu-Mold™ forms to a firm yet pliable digital appliance with no mess or odor. Carefully following the detailed product instructions will produce a perfect digital appliance in just two to three minutes.

Accu-Mold™ has the following characteristics:

- Hypoallergenic, silicone based material.
- Once formed, retains original size and shape.
- Device can be cleansed with soap and water.
- Easily reshaped by grinding or trimming.

Accu-Mold™ is guaranteed for a shelf life of two years.

Medicaments



Clinically proven and popular with Podiatrists, Lumicain® is a colorless hemostatic agent that provides fast and effective control of minor hemorrhage and granulation tissue formation.

9045010 Lumicain®



For the safe and effective treatment of verruca, Salinocaine contains a combination of salicylic acid and an anesthetic for patient comfort. Salinocaine is approved for topical use only.

9045015 Salinocaine

Nail Instruments


Made In England



Nail Nipper, SS, 5 1/2"
• 9085323 (Straight Jaw)



Nail Nipper, SS, 5 1/2"
• 9085324 (Curved Jaw)



Anvil Nail Split Forceps, SS, 5 1/2"
• 9085334

Littauer Cutting Forceps, SS, 4" • 9085393

Littauer Cutting Forceps, SS, 5" • 9085395

Littauer Cutting Forceps, SS, 6" • 9085397



Made In Germany

Collar-Crown 5" • 9085346 (Curved/Sharp)

Collar-Crown 5" • 9085347 (Curved/Blunt)

Curved/Sharp Shown



Collar-Crown 4"
• 9085344 (Curved/Sharp)



Collar-Crown 5" • 9085342 (Straight/Sharp)




Collar-Crown 4"
• 9085345 (Curved/Blunt)



Collar-Crown 4" • 9085340 (Straight/Sharp)

Collar-Crown 4" • 9085341 (Straight/Blunt)

Straight/Sharp Shown



General and Miscellaneous

Needle Holders


Scissors

Retractors


Suction Tubes

Surgical Headband

Instrument Cleaner and Protectant



Needle Holder



Crile-Wood 6" (15.0 cm) TC • 9085153

Tungsten Carbide (TC) Inserts



All Premier needle holders with tungsten carbide inserts, as indicated by a gold handle, are designed for a superior grip, virtually eliminating slippage.



Crile-Wood 6" (15.0 cm) • 9085163



Castroviejo • 9085159

Needle Holder



Mayo-Hegar 6 1/4" (15.6 cm) TC • 9085151



Tungsten Carbide (TC) Inserts



All Premier needle holders with tungsten carbide inserts, as indicated by a gold handle, are designed for a superior grip, virtually eliminating slippage.



Mayo-Hegar 6 1/4" (15.6 cm) • 9085162



Mayo-Hegar 5 1/2" (13.8 cm) • 9085161

Needle Holder



Derf 5" (12.0 cm) TC • 9085155



Derf 5" (12.0 cm) • 9085165



Baumgartner 5" (12.0 cm) • 9085167

Tungsten Carbide (TC) Inserts



All Premier needle holders with tungsten carbide inserts, as indicated by a gold handle, are designed for a superior grip, virtually eliminating slippage.



Baumgartner 5" (12.0 cm) TC • 9085157

Scissors/Nippers

Tissue

Iris - Curved/Delicate 4 1/2" (11.1 cm) • 9085204



Iris - Straight/Delicate 4 1/2" (11.1 cm) • 9085200



Mayo - Straight/Blunt 5 1/2" (13.8 cm) • 9085218



Mayo - Curved/Blunt 5 1/2" (13.8 cm) • 9085220



Scissors/Nippers

Suture



Littauer 4 1/2" (11.1 cm) • 9085114




Littauer 5 1/2" (13.8 cm) • 9085116



Spencer 3 1/2" (8.1 cm) • 9085118



Sawyer Rectal Retractor • 1050105 (Small)
(not actual size)



Sawyer Rectal Retractor • 1050106 (Medium)
(not actual size)

French Scale

	mm	inch
○ 3	1.00	.039
○ 4	1.33	.052
○ 5	1.67	.066
○ 6	2.00	.079
○ 7	2.33	.092
○ 8	2.67	.105
○ 9	3.00	.118
○ 10	3.33	.131
○ 11	3.67	.144
○ 12	4.00	.158

Suction Tubes (not actual size)



(7FR shown)

Suction Tube, Frazier 3FR • 1031203 (with Stylet)

Suction Tube, Frazier 4FR • 1031204 (with Stylet)

Suction Tube, Frazier 5FR • 1031205 (with Stylet)

Suction Tube, Frazier 6FR • 1031206 (with Stylet)

Suction Tube, Frazier 7FR • 1031207 (with Stylet)

Suction Tube, Frazier 8FR • 1031208 (with Stylet)

Suction Tube, Frazier 9FR • 1031209 (with Stylet)

Suction Tube, Frazier 10FR • 1031210 (with Stylet)

Suction Tube, Frazier 11FR • 1031211 (with Stylet)

Suction Tube, Frazier 12FR • 1031212 (with Stylet)



Webster Infusion Cannula • 1050110 (Straight)



Webster Infusion Cannula
• 1050112 (Curved)

Surgical Headband



Jefferson • 1050260

Cleaner

Brite Shield™

Static Soak & Ultrasonic Cleaner

"Most Highly Rated"

Cleans faster Enzymatic action is up to 35% faster than other brands.

Versatile Chemo-mechanically debrides clogged diamond burs, removes alginate and stone from impression trays. Also removes set glass ionomer, carboxylate and zinc phosphate cements. Ideal for removing bio-burden from hand instruments.

Polishes and protects A protective shield coats the instrument through a unique passivation process to restrict rust and corrosion while imparting a lustrous surface to stainless-steel instruments.

The most versatile and economical enzymatic cleaner on the market.

Independent Research Institute



9011100 Case of six - 800gm container with scoop dispenser

Individual units available through your dealer.



Notes

Corporate Headquarters



Our worldwide corporate headquarters, located 10 miles northwest of Philadelphia, PA in Plymouth Meeting, houses marketing/sales for the medical and dental products divisions, customer service, administration and non-manufacturing operations.

Factory



Premier's state-of-the-art factory is located in Philadelphia, PA.

Capabilities include:

- Metal working and finishing
- Injection molding
- Small equipment design, fabrication and servicing
- Product development

If you have a new product idea in the medical or dental field, we encourage you to contact us to discuss possible manufacturing and marketing arrangements.