

Glass Bead Mirror Warmer

Warms Mirrors
Prior to Laryngeal
Examinations to
Prevent Fogging*

- Convenient
- Space-Saving
- UL® Listed



Premier® Medical Products
www.premusa.com

Glass Bead Mirror Warmer

Affording the ENT Clinician improved and efficient laryngeal examinations, while offering the following benefits:

- Rapid warm-up¹ – Ready for use with the day's first patient
- Thermostatically controlled temperature – Turn on in the morning and do nothing until leaving for the evening
- Insulated thermal walls – Housing of the unit remains cool and safe to touch
- Low profile, space-saving design – Less clutter in the examination room
- Large bead well – Less spillage, improved stability
- Hinged lid – Convenient access to bead well, limits exposure to beads when not in use
- UL[®] listed – Meets or surpasses most institutional safety standards



Simply immerse a sterile mirror in the hot beads for three to five seconds prior to examination.

Cautions:

- Allow mirror sufficient time to cool before use
- Soiled or contaminated beads must be changed

1006131 **Mirror Warmer (110V) and Glass Beads**

1006132 **Mirror Warmer³ (220V) and Glass Beads**

1006105 **Replacement Beads**

Designed to fit in your existing cabinet.²

Contact your examination cabinet supplier or surgical supply distributor for details.



Available from:



Premier[®] Medical Products

1710 Romano Drive, Plymouth Meeting, PA 19462 U.S.A.
610-239-6053 • 888 PREM USA (773-6872) • Fax: 610-239-6171
medinfo@premusa.com • www.premusa.com

* UL is a registered trademark of Underwriters Laboratory. / * Not approved for use as a sterilizer. / 1. Inside well reaches a temperature of 250°F after 20 minutes. / 2. Electrical connections should be made only by fully qualified professionals. / 3. The 220V version of the Mirror Warmer is not UL[®] listed.