



Premier Medical Products

# Nitrospray Plus

## Liquid Nitrogen Cryosurgical Instrument

### Effective and efficient treatment of:

- Many benign and premalignant lesions
- Selected malignant lesions\*



Curved trigger provides unobstructed view of treatment site. Trigger can be rotated for left-handed users.

Relief valve precisely regulates operating pressure to deliver a consistent spray and predictable freezing.

Economical luer lock connection. Spray tips (see below) provide versatility and excellent control of nitrogen spray.

Easy-to-hold grooved head fits comfortably in the operator's hand while providing a secure pressure seal.

Stainless construction provides long-term durability and use.

Base holder allows unit to rest securely when not in use.

Spray tip protector (not shown) prevents accidental damage and contact with patient.



### Available in two sizes:

Nitrospray Plus (16oz.) and Nitrospray Plus Lite (10oz.)

All Nitrospray Plus units are backed by a unconditional 30-day and full 3-year warranty.

\* Cryosurgery is not indicated for operative malignant melanoma.