



P.B.S.® Diamond Burs

Nail Reduction & Tissue Debridement

Professional treatment system for hyperkeratotic tissue and hypertrophic nails. The exclusive P.B.S.® bonding system permanently bonds diamonds to the bur surface, providing unsurpassed product life and performance. P.B.S.® burs prove more efficient than steel or carbide instruments and generate less heat for greater patient comfort.

All diamond burs fitted with 3/32'' shaft size.



(To reduce thick, deformed nails.)



(To smooth the rough edges of nails.)



Bud 2026015 (The most popular pattern for nail reduction.)

Short 2026020

(For better control in treating

(For the treatment of incurvated nails or the nail groove.)

Umbrella 2026035

Includes one each:

Bud, Short, Taper, Umbrella

(For the reduction of hyperkeratoses, tylomas and heel callus.)

Diamond Bur System Kit 2026050 A convenient kit of four nail and callus

skin burs plus technique information.



(For working smaller areas.)



(For gross reduction of large nail areas.)



P.B.S. is a registered trademark of the product manufacturer. Abrasive Technology, Inc.



Medicaments



Lumicain®

Clinically proven and popular with Podiatrists, Lumicain® is a colorless hemostatic agent that provides fast and effective control of minor hemorrhage and granulation tissue formation.

9045010 60cc bottle



Salinocaine™

For the safe and effective treatment of verruca, Salinocaine contains a combination of salicyclic acid and an anesthetic for patient comfort. Salinocaine is approved for topical use only.

9045015 40g jar

Uni-Punch®

Seamless Disposable Biopsy Punch

Nine Different Sizes!

Вох	of	25

Size (in mm)	ltem #
1.5	9033509
2	9033502
2.5	9033511
3	9033503
3.5	9033510
4	9033504
5	9033505
6	9033506
8	9033508
Asst.	9033525

Box of 50

Size (in mm)	Item #
1.5	9033519
2	9033512
2.5	9033521
3	9033513
3.5	9033520
4	9033514
5	9033515
6	9033516
8	9033518





