

Seamless Precision



Uni-Punch® Seamless Disposable Biopsy Punch

- ✓ Seamless Blade
Sharpened on both inside and outside
- ✓ Ergonomic Handle Design
*Lightweight, stiff plastic material,
ribbed handle for improved control*
- ✓ Sterile Product
Ethylene oxide sterilization (ETO)
- ✓ Mylar and Tyvek Sterile Peel-Pouch
Individually wrapped
- ✓ Capped Tip for Added Safety
During handling and disposal
- ✓ Packaging Convenience
Boxes of 25 and 50, Assorted package

Available through your authorized dealer.



premiermedicalco.com/promise

premier®

Premier® Medical Products • 888-670-6100 • premiermedicalco.com

©2025 Premier Medical | All Rights Reserved

Uni-Punch®

Seamless Disposable Biopsy Punch



(Assorted - 5 each: 2mm, 3mm, 3.5mm, 4mm, 5mm)

Box of 25

Size (in mm)	Item #
2	9033502
2.5	9033511
3	9033503
3.5	9033510
4	9033504
5	9033505
6	9033506
8	9033508
Asst.	9033525

Box of 50

Size (in mm)	Item #
1.5	9033519
2	9033512
2.5	9033521
3	9033513
3.5	9033520
4	9033514
5	9033515
6	9033516
8	9033518

Nine Different Sizes!
The largest number available!

1.5mm



2mm



2.5mm



3mm



3.5mm



4mm



5mm



6mm



8mm



*Capped tip for
added safety*

premier®



Premier® Medical Products • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premiermedicalco.com • E-mail: info@premiermedicalco.com