

Podiatry Solutions



P.B.S.® Diamond Burs

Nail Reduction & Tissue Debridement

Professional treatment system for hyperkeratotic tissue and hypertrophic nails. The exclusive P.B.S.® bonding system permanently bonds diamonds to the bur surface, providing unsurpassed product life and performance. P.B.S.® burs prove more efficient than steel or carbide instruments and generate less heat for greater patient comfort.

All diamond burs fitted with $\frac{3}{32}$ " shaft size.



Barrel 2026005
(To reduce thick, deformed nails.)



Short 2026020
(For better control in treating smaller digits.)



Umbrella 2026035
(For the reduction of hyperkeratoses, tylomas and heel callus.)



Pear 2026010
(To smooth the rough edges of nails.)



Olive 2026025
(Similar to the PEAR and BUD.)



Small Olive 2026040
(For working smaller areas.)



Bud 2026015
(The most popular pattern for nail reduction.)



Taper 2026030
(For the treatment of incurvated nails or the nail groove.)



Small Cylinder 2026045
(For gross reduction of large nail areas.)

* P.B.S. is a registered trademark of the product manufacturer, Abrasive Technology, LLC.



premiermedicalco.com/promise

premier®

Premier® Medical Products • 888-670-6100 • premiermedicalco.com

©2025 Premier Medical | All Rights Reserved

Medicaments



Lumicain®

Topical Hemostatic Solution

Clinically proven and popular with Podiatrists, Lumicain® is a colorless hemostatic agent that provides fast and effective control of minor hemorrhage and granulation tissue formation.

9045010 60cc bottle



Salinocaine™

Ointment

For the safe and effective treatment of verruca, Salinocaine contains a combination of salicylic acid and an anesthetic for patient comfort. Salinocaine is approved for topical use only.

9045015 40g jar

Uni-Punch®

Seamless Disposable Biopsy Punch

Nine Different Sizes!

Box of 25

Size (in mm)	Item #
1.5	9033509
2	9033502
2.5	9033511
3	9033503
3.5	9033510
4	9033504
5	9033505
6	9033506
8	9033508
Asst.	9033525

Box of 50

Size (in mm)	Item #
2	9033512
2.5	9033521
3	9033513
3.5	9033520
4	9033514
5	9033515
6	9033516
8	9033518



Size 4 shown

Available through your authorized dealer.



Premier® Medical Products • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premiermedicalco.com • E-mail: info@premiermedicalco.com