

ENT



## Glass Bead Mirror Warmer

from Premier® Medical

Warms Mirrors Prior to Laryngeal  
Examinations to Prevent Fogging

**Simply immerse a sterile mirror in the hot beads  
for three to five seconds prior to examination.**

Cautions:

- Allow mirror sufficient time to cool before use
- Soiled or contaminated beads must be changed

*Available through your authorized dealer.*



[premiermedicalco.com/promise](http://premiermedicalco.com/promise)

**premier®**

Premier® Medical Products • 888-670-6100 • [premiermedicalco.com](http://premiermedicalco.com)

©2026 Premier Medical | All Rights Reserved

Affording the ENT Clinician improved and efficient laryngeal examinations, while offering the following benefits:

- ✓ **Rapid warm-up<sup>1</sup>** – Ready for use with the day's first patient
- ✓ **Thermostatically controlled temperature** – Turn on in the morning and do nothing until leaving for the evening
- ✓ **Insulated thermal walls** – Housing of the unit remains cool and safe to touch
- ✓ **Low profile, space-saving design** – Less clutter in the examination room
- ✓ **Large bead well** – Less spillage, improved stability
- ✓ **Hinged lid** – Convenient access to bead well, limits exposure to beads when not in use
- ✓ **UL® listed** – Meets or surpasses most institutional safety standards
- ✓ **Designed to fit in your existing cabinet.<sup>2</sup>** – Contact your examination cabinet supplier or surgical supply distributor for details.



\* UL is a registered trademark of Underwriters Laboratory. / \* Not approved for use as a sterilizer. / 1. Inside well reaches a temperature of 250°F after 20 minutes. 2. Electrical connections should be made only by fully qualified professionals.

1006131 Mirror Warmer (110V) and Glass Beads

1006105 Replacement Beads

Dimensions: 3-1/4"x 5-5/8"x 3-1/4"

Three (3) year warranty



Premier® Medical Products • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA  
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • [premiermedicalco.com](http://premiermedicalco.com) • [info@premiermedicalco.com](mailto:info@premiermedicalco.com)