

DERMATOLOGY

Nitrospray Plus

Liquid Nitrogen Cryosurgery Instrument

Trigger Operated, Closed System, Liquid Nitrogen Cryosurgical Unit, while offering the following benefits:

- ✓ Grooved Head – *Easy to use, comfortable to hold.*
- ✓ Curved Trigger – *Unobstructed view of treatment site.*
- ✓ Notched Head Grip – *Easy to remove/install cap.*
- ✓ Three (3) year warranty

Effective and efficient treatment of:

- Many benign and premalignant lesions
- Selected malignant lesions*

Static holding time: 6 – 8 hours

Intermittent use: 3 – 6 hours of patient treatment time

Available in two sizes:

Nitrospray Plus (16oz.) and Nitrospray Plus Lite (10oz.)

All Nitrospray Plus units are backed by a unconditional 30-day and full 3-year warranty.

1006060 Nitrospray Plus (16oz fill capacity)

1006065 Nitrospray Plus Lite (10oz fill capacity)

Includes: 5 spray tips, Tip protector, Base holder

*Please see our website for full listing of Nitrospray Plus accessories.

*Cryosurgery is not indicated for operative malignant melanoma.



premiermedicalco.com/promise

premier

Premier® Medical Products • 888-670-6100 • premiermedicalco.com

0323039 Rev2